

CerviGrip™

Uterine Manipulator

Secure Cervical Fixation with Adjustable Uterine Tip



General Laparoscopy

Adjustable uterine tip to fit all uterine lengths ✓

Built-in tenaculum holder ✓

Tight cervical seal for hydrotubation ✓

Syringe port for hydrotubation ✓

TLH / LAVH

Adjustable tip removes need for various tips ✓

Secure cervical fixation ✓

Predictable ureteric and uterine distances with ColpoWave ✓

Facilitates uterine removal ✓



SurgiTools.com